



Recognising right and left hearing aid

Your hearing aids are individually tuned, therefore it is important you wear the correct hearing aid in each ear. Do not swap them. Please pay attention to this when cleaning, storing and inserting your hearing aids.

The red marker indicates the right hearing aid and the blue marker the left hearing aid.

Depending on the model of hearing aid, your marker will be in one of the following places.

If you cannot see any coloured markers please contact your local audiology department.



Audiologist contact details:

The right choice for the NHS

UK 08.06 A 01.2020

danalogic GN