

# Recognising right and left hearing aids


Your hearing aids are programmed for each ear. It is important that you wear the correct hearing aid in each ear and do not swap them. Your hearing aids have coloured markers to help you know which ear they are programmed for.

The red marker indicates the right hearing aid and the blue marker indicates the left hearing aid.


The images below show you where you can find your coloured marker dependent on the model of your hearing aid.

If you cannot see any coloured markers please contact your local Audiology department.


For 62/77/88 models



For 67 models



For 98 models



For further support, scan the QR code or visit our patient support portal at: [danalogic.co.uk/patient-portal](https://danalogic.co.uk/patient-portal)