


Recognising right and left hearing aids for danalogic Extend



Your hearing aids are programmed for each ear. It is important that you wear the correct hearing aid in each ear and do not swap them. Your hearing aids have coloured markers to help you know which ear they are programmed for.

The red marker indicates the right hearing aid and the blue marker indicates the left hearing aid.

The image on this page shows you where you can find your coloured marker.

If you cannot see any coloured markers please contact your local Audiology department.



For further support, scan the QR code or visit our patient support portal at: danalogic.co.uk/patient-portal

The right choice for the NHS

danalogic GN